





KA 200 N - cardioid



No. 82 - 3020


KA 200 N cardioid

recommended preamp	MBP 603 (A)	MBP 680	MBP 648 (A)
frequency response Hz	40 - 20.000	40 - 20.000	40 - 20.000
sensitivity	14 mV/Pa	7 mV/Pa	7 mV/Pa
signal to noise ratio CCIR	69 dB	63 dB	68 dB
equivalent SPL rated at CCIR	25 dB	31 dB	26 dB
signal to noise ratio DIN / IEC	82 dB	76 dB	81 dB
equivalent SPL rated at DIN / IEC	12 dB	18 dB	13 dB
max. spl	130 dB	130 dB	130 dB
size	ø 21 x 20 mm		



- cardioid pattern
- highly versatile capsule
- maximum rejection on rear-incident sound
- X/Y recordings or Mid-Side configuration
- frequency response 40 - 20.000 Hz
- neutral record space ± 135°


KA 400 N - cardioid



No. 82 - 3040

KA 400 N cardioid

recommended preamp	MBP 603 (A)	MBP 680	MBP 648 (A)
frequency response Hz	50 - 20.000	50 - 20.000	50 - 20.000
sensitivity	15 mV/Pa	7 mV/Pa	7 mV/Pa
signal to noise ratio CCIR	70 dB	63 dB	68 dB
equivalent SPL rated at CCIR	24 dB	31 dB	26 dB
signal to noise ratio DIN / IEC	83 dB	76 dB	81 dB
equivalent SPL rated at DIN / IEC	11 dB	18 dB	13 dB
max. SPL	131 dB	131 dB	131 dB
size	ø 21 x 20 mm		



- cardioid speech optimized
- 12 dB low cut at 50 Hz
- frequency response 50 - 20.000 Hz